

PATCH TESTING : PATIENT TEST RESULTS

Patient Name: _____ ID#: _____

Panel Placement: _____ (date/time) Panel Removal: _____ (date/time)

Read Date #1: _____ Read Date #2: _____ Read Date #3: _____ Read Date #4: _____

PG-1000 Plastics & Glues Series				Readings			
#	Article#	Veh/Conc	Description	#1	#2	#3	#4
1	H-007	1.0 pet	Hydroquinone				
2	D-007	5.0 pet	Dibutyl phthalate				
3	P-011	1.0 pet	Phenyl salicylate(Salol)				
4	D-018	2.0 pet	Diocetyl phtalate(DEHP,DOP) / (Dimethylhexylphthalate)				
5	D-006	2.0 pet	BHT / (2,6-Di-tert-butyl-4-cresol)				
6	H-016	1.0 pet	Drometrizole / (2(2-Hydroxy-5-methyl-phenyl)benzotriazol)				
7	B-007	1.0 pet	Benzoylperoxide				
8	B-030B	0.25 pet	4-tert-Butylcatechol (PTBC)				
9	A-018	1.0 pet	Azodiisobutyrodinitrile				
10	B-013	1.0 pet	Bisphenol A/(4,4'-Isopropylidene-diphenol)				
11	T-015	5.0 pet	Tricresyl phosphate				
12	P-005	1.0 pet	Phenol formaldehyde resin (P-F-R-2)				
13	B-024	1.0 pet	4-tert-Butylphenolformaldehyde resin				
14	T-022	5.0 pet	Triphenyl phosphate				
15	T-010	10.0 pet	Toluenesulfonamide formaldehyde resin				
16	R-002	1.0 pet	Resorcinol monobenzoate				
17	P-007	2.0 pet	2-Phenylindole				
18	B-022	2.0 pet	2-tert-Butyl-4-methoxyphenol (BHA)				
19	A-002	10.0 pet	Hydroabietyl alcohol / (Abitol)				
20	B-023	1.0 pet	4-tert-Butylphenol				
21	M-015	1.0 pet	2-Monomethylol phenol/(Saligenin)				
22	D-025	1.0 pet	Diphenylthiourea				
23	O-004	0.1 pet	2-n-Octyl-4-isothiazolin-3-one				
24	C-027	1.0 pet	Cyclohexanone resin				
25	T-028	0.5 pet	Triglycidyl isocyanurate				

Notes
